

Electrical assistant for your manual wheelchair



Images for illustration purposes only. Certain features may be optional.

## **Technical information V-Drive**

seatwidth	max.	max.	action radius	power	battery	kg weight	users weight
390 - 480 500 - 600	6 km/h	8°	ca. 15 km	150 W 200 W	2x12V/12Ah 2x12V/12Ah (H.D. version)	21 kg 1 <i>7</i> ,4 kg	135 kg* 200 kg*
(H.D. version)				(H.D. version)	(H.D. version)		(H.D. version)

Measurements in mm, weight in kg - Deviations in sizes and weight (1,5 cm and 1,5 kg) are possible. \* wheelchair included.

## **Comfort**

The V-Drive system simplifies your movements with a manual wheel-chair. The V-Drive takes over all the work from your pushing assistant. Taking a steep hill or inclinations,



bumps in the roads or just a small journey outside... all this is just ideal for the V-Drive.

A userfriendly controller ensures optimal quality and comfort for daily usage of your wheelchair in combination with the V-Drive.

We are offering two versions of the V-Drive. A standardversion (for people up to 135 kg) and a Heavy-Duty version (for people up to 200 kg).

## Installation

A very simple installation is available for left or rightside mounting of the V-Drive. Disassembly of your V-Drive for transport inside your car can be done in a minute. Disconnect the battery pack and the driving system and you are ready to go.

If you prefer not to use your V-Drive along the way, it is sufficient to lift up the driving system, hang it to on to one of the pushing handles and you have regained manual control over your wheelchair.

The **NEW** height adjustment allows the motor to be positioned horizontally.



## N.V. Vermeiren N.V.

Vermeirenplein 1-15 B-2920 Kalmthout Belgium

Tel.: +32 (0) 3 620 20 20 Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be

For more information, contact your local dealer



: N.V. Vermeiren N.V., Vermeirenplein 1/15 - 2920 Kalmthout - Belgium -10/2016